

Manual

qPCR-HS Mix SYBR®

High specificity ready-to-use mix for real-time Hot Start PCR with SYBR® Green. Mixture contains monoclonal antibody blocked Taq DNA polymerase (RUN-HS).

catalog #	size
2008HS-100	200 reactions in 25 µl
2008HS-1000	2000 reactions in 25 µl

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.
The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

SYBR® is a registered trademark of Molecular Probes Inc.



Description

qPCR-HS Mix SYBR® is optimized for high specificity ready to use real-time Hot Start PCR mixture with SYBR® Green dye. Mixture contains all components required for qPCR except DNA template and primers. Activation of the monoclonal antibody blocked RUN-HS polymerase occurs during initial denaturation in PCR.

Contents

	2008HS-100		2008HS-1000		storage
	quantity	cat #	quantity	cat #	
2x qPCR-HS Mix SYBR® (qPCR-HS Mix SG)	2 x 1.25 ml	K-28-125A	20 x 1.25 ml	K-28-125A	-20 °C
ultrapure water	2 x 1.5 ml	K-WUP-15A	20 x 1.5 ml	K-WUP-15A	-20 °C

Notes

- Before use, it is necessary to completely thaw and thoroughly mix the kit components by gently inverting the tube.
- Up to 7x repeated freeze-thaw cycles do not influence the activity of this product.

ROX reference dye

Some PCR instruments perform fluorescence signal correction and it is recommended to use ROX reference dye for signal normalization. Please follow manufacturer's instructions regarding addition of ROX reference dye and its concentration.

Example qPCR protocol

1. Add to the PCR tubes:

component	volume	final concentration
	25 µl	
2x qPCR-HS Mix SYBR®	12.5 µl	1X
primer 1 (10 µM)*	0.5 µl	0.2 µM
primer 2 (10 µM)*	0.5 µl	0.2 µM
DNA template	1-5 µl	< 250 ng/reakcja
ultrapure water	up to 25 µl	

*For optimization, a primer titration should be performed from 0,2 µM do 1 µM final concentration.

2. Gently mix the samples and briefly centrifuge.

3. Place the tubes in the thermocycler and start the PCR programme.
An example amplification profile:

reaction step	temperature	time	number of cycles
enzyme activation	95 °C	5 min	1
denaturation	95 °C	15 s	40
annealing*	50-68 °C	30 s	
extension**	72 °C	30 s	

*Annealing temperature depends on primer sequence and the composition of the reaction mixture.

**Time of extension depends on the length of the amplicon.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone +48 883 323 761, +48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2025-1

