

Manual

PCR Kit 5

Complete kit for PCR including *Taq* DNA polymerase and reaction buffers. Concentration 5 U/ μ l.

catalog #	size	concentration
1205-200	200 U	5 U/ μ l
1205-1000	1000 U	5 U/ μ l

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

Advantages

- Complete kit recommended for standard PCR reaction.
- This kit contains the most popular used thermostable *Taq* DNA polymerase.

Description

Taq DNA polymerase is thermophilic DNA polymerase purified from *E.coli* stream carrying a plasmid with a cloned gene encoding a DNA polymerase from *Thermus aquaticus*.

Enzyme catalysis incorporation of deoxynucleotides to 3' end of dsDNA at temp. 70-80 °C and presence of Mg^{2+} ions.

Taq DNA polymerase lacks 3'-5' exonuclease activity (proofreading), but possesses weak 5'-3' exonuclease activity.

KU buffer increases the specificity of the PCR reaction for DNA templates with secondary structures and GC pairs.

Using KU buffer it's necessary to prepare a control reaction without KU buffer.

"I" reaction buffer contains Mg^{2+} ions at a concentration ensuring satisfactory results in most experimental systems.

Optimization of the concentration of Mg^{2+} ions in the reaction is the possibility of using an "III" reaction buffer (without Mg^{2+} ions) and adding an appropriate amount of Mg^{2+} ions to it in the form of $MgCl_2$ included in the kit.

Contents

	1205-200	1205-1000	storage
RUN polymerase	200 U	1000 U	-20 °C
dNTP Mix (10 mM)	200 µl	4 x 200 µl	-20 °C
10x "I" buffer (with Mg^{2+} ions) 100 mM KCl, 100 mM $(NH_4)_2SO_4$, 200 mM Tris-HCl, pH 8.5, 15 mM $MgSO_4$, 1% Triton X-100	1.5 ml	4 x 1.5 ml	-20 °C
10x "III" buffer (without Mg^{2+} ions) 100 mM KCl, 100 mM $(NH_4)_2SO_4$, 200 mM Tris-HCl, pH 8.5, 1% Triton X-100	1.5 ml	4 x 1.5 ml	-20 °C
5x KU buffer no DMSO, no toxic reagents	2 ml	4 x 2 ml	-20 °C
6x loading buffer	1 ml	1 ml	-20 °C
$MgCl_2$ (25 mM)	1.5 ml	2 x 1.5 ml	-20 °C
ultrapure water	5 ml	4 x 5 ml	-20 °C

Example PCR protocol

- 1. Thaw all components on ice, gently mix by inverting the tube and briefly centrifuge. Place the tube on ice.
- 2. Place PCR tube on ice and add:

component	PCR reaction volume
	25 µl
10x “I” buffer or “III” buffer	2.5 µl
dNTP Mix (10 mM)	200-250 µM (0.5-0.6 µl)
primer 1	0.1-0.5 µM
primer 2	0.1-0.5 µM
RUN polymerase	1-2 U
DNA template	10 pg - 1 µg
5x KU buffer (option)	2.5-5 µl
MgCl ₂ (option)	depending on the needs
6x loading buffer (option)	depending on the needs
ultrapure water	up to 25 µl

- 3. Gently mix the sample and briefly centrifuge.
- 4. Place the tube in the thermocycler and start the PCR programme.

An example amplification profile for products up to 1000 bp:

step	temperature	time
initial denaturation	95 °C	1-3 min
25-45 cycles	94 °C	30-60 s
	50-68 °C	30-60 s
	72 °C	1 min
final incubation	72 °C	5-10 min

- 5. PCR product store in a refrigerator or freezer until later use.



A&A BIOTECHNOLOGY

innovating life science

A&A Biotechnology, ul. Strzelca 40, 880-299 Gdańsk, Poland
phone: +48 883 323 761, +48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2024-1

