

## Manual

# Fast DNA Plant Screen PCR

Kit for rapid isolation of genomic DNA from plant material, to be used in PCR.

catalog #	size
050-192P	192 isolations (2 x 96)

For research use only.

### Guarantee

A&A Biotechnology provides a guarantee on this product.

The company does not guarantee the correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

# Table of Contents

<b>Contents</b>	<b>3</b>
<b>Additional equipment and reagents</b>	<b>3</b>
Necessary	3
Optional	3
<b>Isolation protocol</b>	<b>4</b>
<b>Example PCR protocol</b>	<b>5</b>

# Contents

component	quantity	storage
A buffer	5 x 4 ml	-20 °C
B buffer	5 x 4 ml	-20 °C
PCR Mix Plus HGC	1 pcs	-20 °C

## Additional equipment and reagents

### Necessary

- 1,5 ml Eppendorf tubes
- Heatblock or incubator set to 95 °C

### Optional

- PCR tubes or 96-well PCR plates

## Isolation protocol

1. Cut the plant tissue in small fragments (size not exceeding 2-3 mm) and transfer to 1,5 ml Eppendorf tubes (not included).
2. Add **50-100 µl** of **A** buffer. Mix by pipetting. Sample must be completely submerged in A buffer.
3. Incubate the sample for **10 min** at **95 °C**.
4. Cool down the samples to **room temp**.
5. Add an equal volume of **B** buffer (about **50-100 µl**).
6. Mix the samples for **5 s** by pipetting or vortexing.
7. Store the samples with fragments of plant tissue up to **3 months** at **4-8 °C**.  
  
**Do not freeze!**  
Freezing may cause DNA degradation.

We recommend using isolated DNA to maximum 10% of the final volume of PCR sample (e.g. 5 µl of isolated DNA to 50 µl of final volume of PCR mix reaction).

## Example PCR protocol

We recommend using isolated DNA to maximum 10% of the final volume of PCR sample (e.g. 5 µl of isolated DNA to 50 µl of final volume of PCR mix reaction).

- Before use all solutions should be thawed thoroughly on ice, gently mixed by inverting the tube and briefly centrifuged.
- Up to 7x repeated freeze-thaw cycles do not influence the activity of this product.

1. Thaw the **PCR Mix Plus HGC** and **ultrapure water** on ice, gently mix by inverting the tubes and briefly centrifuge. Place the tubes on ice again.
2. Place PCR tubes on ice and add:

component	PCR reaction volume	
	25 µl	50 µl
<b>PCR Mix Plus HGC</b>	12.5 µl	25 µl
<b>Starter 1</b>	0.1-1 µM	0.1-1 µM
<b>Starter 2</b>	0.1-1 µM	0.1-1 µM
<b>DNA template</b>	10 pg-1 µg	10 pg-1 µg
<b>Ultrapure water</b>	up to 25 µl	up to 50 µl

3. Gently mix the samples and briefly centrifuge. If necessary, overlay the samples with mineral oil. (it's recommended for thermocyclers without a heated lid).
4. Place the tubes in the thermocycler and start the PCR programme.

An example amplification profile for products up to 500 bp:

step	temperature	time
<b>Initial denaturation</b>	95 °C	2-3 min
<b>25 - 45 cycles</b>	95 °C	15-30 s
	50 - 68 °C	30-60 s
	72 °C	15-60 s

5. Load the post-PCR samples directly on an agarose gel for electrophoresis.







**A&A BIOTECHNOLOGY**  
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland  
phone +48 883 323 761, +48 600 776 268  
info@aabiotech.com, www.aabiotech.com

version 2023-1

