

Manual

WALK DNA polymerase

Pwo DNA polymerase with reaction buffer. Concentration 1 U/ μ l.

catalog #	size
1002-200	200 U
1002-1000	1000 U

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

Advantages

- Highest fidelity thermostable *Pwo* DNA polymerase (*Pyrococcus woesei*)
- Proofreading
- Highest purity enzyme
- Recommended for cloning - generates blunt-ended DNA

Description

WALK DNA polymerase is recombinant thermophilic DNA polymerase purified from *E.coli* stream carrying a plasmid with a cloned gene encoding a DNA polymerase from *Pyrococcus woesei*.
Enzyme catalysis incorporation of deoxynucleotides to 3' end of dsDNA at temperature 70-80 °C and presence of Mg²⁺ ions. Enzyme generates blunt-ended ds DNA fragments.
Pwo DNA polymerase possesses 3'-5' exonuclease activity (proofreading) responsible for high fidelity of the enzyme, but lacks 5'-3' exonuclease activity.

Contents

	1002-200	1002-1000	storage
WALK polymerase	200 U (1 U/μl)	1000 U (1 U/μl)	-20 °C
storage buffer: 10 mM KCl, 20 mM Tris-HCl pH 8.7, 0,1 mM EDTA, stabilizers, 50% glicerol (v/v).			
WALK reaction buffer	1 x 1.5 ml	4 x 1.5 ml	-20 °C
10x PCR reaction buffer: 100 mM KCl, 100 mM (NH ₄) ₂ SO ₄ , 200 mM Tris-HCl, pH 8.5, 20 mM MgSO ₄ , 1% Triton X-100.			

Notes

- Before using, thoroughly thaw and gently mix by inverting the tubes.

Example PCR protocol

1.
- Thaw all components on ice, gently mix by inverting the tubes and briefly centrifuge. Place the tubes on ice again.
2.
- Place PCR tubes on ice and add:

component	PCR reaction volume
	50 µl
WALK reaction buffer	5 µl
dNTP Mix (10 mM)	200-250 µM (1-1.25 µl)
Starter 1	0,1-0,5 µM
Starter 2	0,1-0,5 µM
WALK polymerase	1-2 U
DNA template	10 pg - 1 µg
Sterile water	up to 50 µl

3.
- Gently mix the samples and briefly centrifuge.
4.
- Place the tubes in the thermocycler and start the PCR programme.

An example amplification profile for products up to 1000 bp:

step	temperature	time
Initial denaturation	94 °C	1-5 min
25-45 cycles	94 °C	30-60 s
	50-68 °C	30-60 s
	72 °C	2 min
Final incubation	72 °C	5-10 min

For longer fragments we recommend doubling the elongation time in relation to the time used for *Taq* DNA polymerase (RUN DNA polymerase).

5.
- PCR products store in a refrigerator or freezer until later use.



A&A BIOTECHNOLOGY

innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone: +48 883 323 761, +48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2023-2

