

Manual

RT HS-PCR Mix SYBR®

Ready-to-use mix for real-time hot-start PCR with SYBR® Green.
2x concentrated.

catalog #	size
2017-100HS	200 reactions in 25 µl
2017-1000HS	2000 reactions in 25 µl

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

SYBR® is a registered trademark of Molecular Probes Inc.

Description

RT HS-PCR Mix SYBR® are optimized ready to use real-time hot-start PCR mixtures with SYBR® Green. Mixture contains all components required for qPCR except DNA template and primers.

This mix is recommended especially for difficult PCRs and multiplex PCR. StartWarm technology allows for activation of *Taq* DNA polymerase in the first steps of PCR. Long time of preliminary denaturation step is not required.

The premix formulation saves time and reduces contamination due to a reduced number of pipetting steps required for PCR set up. The mix is optimized for efficient and reproducible reaction.

Contents

	2017-100HS	2017-1000HS	storage
RT HS-PCR Mix SYBR®	2 x 1.25 ml	20 x 1.25 ml	-20 °C, in darkness
ultrapure water	2 x 1.5 ml	20 x 1.5 ml	-20 °C

RT HS-PCR Mix SYBR® A/B/C composition

component	amount
<i>Taq</i> DNA polymerase	0.1 U/µl
MgCl ₂	4 mM
dNTPs	0.5 mM of each dNTP
2x reaction buffer with SYBR® Green	

Notes

- Before use all solutions should be thawed thoroughly on ice, gently mixed by inverting the tube and briefly centrifuged.
- Up to 7x repeated freeze-thaw cycles do not influence the activity of this product.

Example PCR protocol

1.
- Thaw **all components of the kit** on ice, gently mix by inverting the tubes and briefly centrifuge. Place the tubes on ice again.
2.
- Place PCR tubes on ice and add:

component	PCR reaction volume		
	10 µl	25 µl	50 µl
RT HS-PCR Mix SYBR®	5 µl	12.5 µl	25 µl
primer 1**	0.1-1 µM*	0.1-1 µM*	0.1-1 µM*
primer 2**	0.1-1 µM*	0.1-1 µM*	0.1-1 µM*
DNA, cDNA template	10 pg-1 µg	10 pg-1 µg	10 pg-1 µg
ultrapure water	up to 10 µl	up to 25 µl	up to 50 µl

* recommended for standard qPCR
** final concentration in reaction mixture

3.
- Gently vortex the samples and briefly centrifuge to collect all droplets remaining on the tube walls and caps to the bottom of the tube.
4.
- Place the tubes in the thermocycler and start the PCR programme.

An example amplification profile:

step	temperature	time
initial denaturation	95 °C	3-5 min
25-45 cycles	95 °C	15-30 s
	50-68 °C	30-60 s
	72 °C	15-60 s*

* depending on the length of PCR products
PCR product melting analysis is recommended.

Recommended ROX mixture

HiROX (0.6-1 µl per 50 µl of total reaction volume): Applied Biosystems: 7000, 7300, 7700, 7900HT Fast, StepOne, StepOnePlus.

LowROX (0.6-1 µl per 50 µl of total reaction volume): Applied Biosystems: 7500, Stratagene: Mx3000P, Mx3005P, Mx4000P.



A&A BIOTECHNOLOGY

innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone +48 883 323 761, +48 600 776 268
info@aabiotech.com, www.aabiotech.com

