

## Manual

# MagnifiQ™ 1 Pathogen instant kit

Kit for automated, magnetic isolation of DNA and RNA of pathogens in the strip format. Contains ready-to-use, reagent-filled stripes and all necessary consumables.

catalog #	size	compatible devices *
607A-1V-32	32 isolations	Auto-Pure Mini, Auto-Pure S32
607A-1V-160	160 isolations	Auto-Pure Mini, Auto-Pure S32

#### \* Compatible devices

The kit has been tested with specific Allsheng brand isolation devices. This does not preclude it from working with other devices. If your device is not listed, please contact us at [info@aabiotech.com](mailto:info@aabiotech.com).

For research use only.

#### Guarantee

A&A Biotechnology provides a guarantee on this product.

The company does not guarantee the correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components.

# Table of Contents

<b>Advantages</b>	<b>3</b>
<b>Sample types</b>	<b>3</b>
<b>Specification</b>	<b>3</b>
<b>Description</b>	<b>3</b>
<b>Contents</b>	<b>3</b>
<b>Additional equipment and reagents</b>	<b>4</b>
Necessary	4
Optional	4
<b>Material preparation</b>	<b>4</b>
Animal Tissue	4
Blood (fresh or frozen, plasma, serum), Body Fluids	4
Swabs with transport medium	5
Dry swabs	5
<b>Protocol</b>	<b>6</b>
Protocol files	6
Extraction Protocol	7
<b>Safety information</b>	<b>9</b>

## Advantages

- Requires only a few minutes of manual work while adding the samples. The rest of the procedure is carried out in an automatic extraction device.
- It does not require the initial preparation of buffers. Just add prepared samples to the plate and after about 30 min you get extracted material.
- Universal extraction of both DNA and RNA.

## Sample types

sample type	sample size
<a href="#">Animal tissue</a>	up to 20 mg
<a href="#">Blood (fresh or frozen, plasma, serum), body fluids</a>	up to 200 µl
<a href="#">Swab</a>	1 pcs

## Specification

protocol time	~ 30 min
elution volume	100 µl
elution solution	Tris buffer (10 mM, pH 8.5)
binding capacity	30 µg DNA/RNA
downstream applications	qPCR, RT-qPCR, sequencing

## Description

**MagnifiQ™ 1 Pathogen instant kit** is designed for RNA and DNA isolation from viruses and Gram(-) bacteria. The isolated material is perfect for further analyzes and tests by qPCR and RT-PCR methods and for sequencing.

The **MagnifiQ™** product series is based on the automated isolation of nucleic acids with use of magnetic beads. This method significantly shortens working time and reduces risk of mistakes in comparison to manual methods.

## Contents

component	607A-1V-32		607A-1V-160		storage
	quantity	cat #	quantity	cat #	
XS-G - extraction strip	4 x 8 pcs	K-S1V22XG	20 x 8 pcs	K-S1V22XG	15–25 °C
LTE 2X buffer	8 ml	K-LTE2X-8	35 ml	K-LTE2X-35	15–25 °C
Proteinase K	1,1 ml	K-PRK-11A	4 x 1,1 ml	K-PRK-11A	2–8 °C*
tip comb 8	16 pcs	K-C8U-16	2 x 40 pcs	K-C8U-40	15–25 °C

\* Proteinase K can be stored at 15–25 °C for up to 12 months.

# Additional equipment and reagents

## Necessary

- pipette
- pipette tips
- 1.5 ml Eppendorf tubes with lock (sample lysis)

## Optional

- vortex
- sterile water, Tris buffer, PBS buffer

## Material preparation

### Animal tissue

1. Homogenize up to **20 mg** of animal tissue in the PBS or Tris buffer.
2. Centrifuge the sample for **1 min** at **500 RPM**.
3. Transfer **200 µl** of the supernatant to 1.5 ml Eppendorf tube with lock.
4. Add **200 µl** of **LTE 2X** buffer and **20 µl** of **Proteinase K**.
5. Vortex the sample for **10 s**.
6. Follow point 1. of the [Extraction Protocol](#).

### Blood (fresh or frozen, plasma, serum), body fluids

1. Transfer **200 µl** of the sample to 1.5 ml Eppendorf tube with lock.
2. Add **200 µl** of **LTE 2X** and **20 µl** of **Proteinase K**.
3. Vortex the sample for **10 s** and incubate for **10 min** at **50 °C** with shaking.  
**Note.** If automatic shaking is not available, mix the samples by inverting the tubes a few times.
4. Centrifuge for **20 s** at **10 000 RPM**.

5. **Attention.** In the Extraction Protocol, use the supernatant as the sample.  
Follow point 1. of the [Extraction Protocol](#).

## Swabs with transport medium

No additional material preparation is required.




- 1 Follow point 1. of the [Extraction Protocol](#).

## Dry swabs

1. Break or cut off part of the swab with the collected sample and place it in 1.5 ml Eppendorf tube with lock.  
**Note.** The portion of the swab with the collected sample should fit completely into the tube.
2. Add **500 µl** of sterile water, **Tris** buffer or **PBS** buffer.  
**Note.** Part of the swab with the sample should be completely immersed in the buffer.
3. Leave at room temperature for **10 min**.
4. Vortex for **10 s**.
5. Transfer **200 µl** of supernatant to new 1.5 ml Eppendorf tube with lock.
6. Add **200 µl** of **LTE 2X** buffer and **20 µl** of **Proteinase K**.
7. Follow point 1. of the [Extraction Protocol](#).


# Protocol

## Protocol files

device	protocol name	protocol file	installation
Auto-Pure Mini	MQ-UNI-MI	<a href="http://aabiotech.com/protocols/magnifiq/MI/MQ-UNI-MI.txt">aabiotech.com/protocols/magnifiq/MI/MQ-UNI-MI.txt</a>	<ol style="list-style-type: none"> <li>1. Create folder "items" on a USB drive and copy the protocol file to it.</li> <li>2. Insert the USB drive into a USB slot in the device.</li> <li>3. On a device screen, go to Settings &gt; System &gt; Transfer &gt; Import.</li> <li>4. Select the protocol and tap "Import".</li> </ol>
Auto-Pure Mini (QR code)	MQ-UNI-MI		<ol style="list-style-type: none"> <li>1. On a device screen, go to Run &gt;  &gt; </li> <li>2. Scan the QR code with the device's scanner.</li> </ol>
Auto-Pure S32	MQ_UNI_S32	<a href="http://aabiotech.com/protocols/magnifiq/S32/MQ_UNI_S32.txt">aabiotech.com/protocols/magnifiq/S32/MQ_UNI_S32.txt</a>	<ol style="list-style-type: none"> <li>1. Create folder "im_export_protocols" on a USB drive and copy the protocol file to it.</li> <li>2. Insert the USB drive into a USB slot in the device.</li> <li>3. On a device screen, go to Protocols &gt; Import.</li> <li>4. Select the protocol and tap "Import".</li> </ol>

## Extraction Protocol

1. Place **XS-G** stripes in the rack(s).



2. Remove the foil from the **XS-G** stripes.  
**Note:** It is important to do this carefully so as not to mix buffers from different wells.  
**Note:** If possible, briefly spin before removing the foil.
3. Add **400 µl** of sample to the well **1** (first from the left) on the **XS-G** strip.  
**Note:** The wells are numbered on the side of the strip.
4. Place the racks(s) in the extraction device.
5. Place the appropriate number of **tip combs 8** in the extraction device.
6. Run the protocol on your device.

7. After the run is over, first remove the combs and then the racks from the extraction device and transfer the purified DNA / RNA located in well **6** (first from the right) on the **XS-G** strip to a sterile tube.  
  
**Note.** For longer storage of extracted material, transfer it from the plate to appropriate tubes and store at 4 °C for DNA or -70 °C for RNA.



# Safety information



**DANGER**

## Proteinase K


H315 Causes skin irritation.  
 H319 Causes serious eye irritation.  
 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
 H335 May cause respiratory irritation.  
 P261 Avoid breathing dust.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P342+P311 If experiencing respiratory symptoms call a Poison Center or doctor/physician.



**WARNING**

## LTE 2X

H302 Harmful if swallowed.  
 H315 Causes skin irritation.  
 H319 Causes serious eye irritation.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



**DANGER**

## XS-G extraction strip

H225 Highly flammable liquid and vapor.  
 H302+H312+H332 Harmful if swallowed, in contact with skin or if inhaled.  
 H314 Causes severe skin burns and eye damage.  
 H412 Harmful to aquatic life with long lasting effects.  
 P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
 P273 Avoid release to the environment.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.  
 P301+P312+P330 If swallowed: Call a poison center/doctor/ if you feel unwell.  
 P303+P361+P353 If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.  
 P304+P340 If inhaled: Remove person to fresh air and keep comfortable for breathing.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.







**A&A BIOTECHNOLOGY**  
innovating life science

A&A Biotechnology, Strzelca 40, 80-299 Gdańsk  
tel. 883 323 761, 600 776 268  
info@aabiotech.com, www.aabiotech.com

version 2025-1

