

## Manual

# 6x loading buffer for DNA electrophoresis

6x concentrated DNA sample loading buffer for electrophoresis. Dyes: bromophenol blue, xylene cyanole.

cat #	size
5000-5	5 x 1 ml

For research use only.

### Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components



# Contents

## Loading buffer for DNA electrophoresis

- 20% ficoll
- 0.02% bromphenol blue
- 0.02% xylene cyanole FF

## Storage

Store at 2–8 °C

## Application

Use 1 µl of loading buffer per 5 µl of sample.



# A&A BIOTECHNOLOGY

innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland  
phone +48 883 323 761, +48 600 776 268  
info@aabiotech.com, www.aabiotech.com

version 2023-1

