



ISO 13485

LL-C (Certification)




MagnifiQ Blood Direct

work smarter, not harder



A&A BIOTECHNOLOGY
innovating life science

Simple, fast, and automated DNA isolation from blood



Add blood sample

Start automated purification

Obtain high-quality DNA

Advantages

- Simple and fast, automated isolation procedure
- No sample pre-treatment required
- Complete sample isolation performed directly within the automated system
- Repeatable results
- High purity nucleic acids suitable for downstream applications

Specification

protocol time	~ 45 min
sample type	blood fresh or frozen
sample size	up to 200 µl
elution volume	100 µl
elution solution	Tris buffer (10 mM, pH 8.5)
binding capacity	30 µg DNA/RNA

Compatible devices


Auto-Pure Mini • Auto-Pure S32 • Auto-Pure 32A

Don't see your device on the list? Contact us orders@aabiot.com – we're happy to help!

Ordering information

MagnifiQ™ 1 Blood Direct instant kit - strip format

MagnifiQ™ 16 Blood Direct instant kit - 16 samples per plate format

Find out more  www.aabiot.com! Contact us - let's create the perfect solution for you.



Produced in the EU – Gdańsk, Poland

 www.aabiot.com
 orders@aabiot.com
 +48 600 776 268
 Strzelca 40
80-299 Gdańsk